



## VITROS<sup>®</sup> XT Solutions

Consistently:  
Fast. Accurate. Reliable.

Because **Every**  
Test is a **Life**<sup>™</sup>



Every day, lives depend on your lab's results.

**You're under constant pressure to deliver the best patient care.**


**How often:**

- Do you question a critical result?
  - Are you forced to repeat a test?
  - Is your system not available?
  - Do you meet your stat turnaround time?
- 

**Do you trust your current chemistry/immunoassay system to consistently deliver?**

# Meet VITROS® XT Systems: Patient care, delivered consistently fastest and more reliably accurate.

## Efficiency




Solutions that meet your operational and cost efficiency objectives

## Confidence



Confidence in the accuracy and treatment impact of results you deliver

## Reliability



A system you can rely on to be ready to deliver accurate results when you need them

## Ownership



Fully supported and worry free



**The VITROS® XT 3400 Chemistry System**

Up to 1,130 tests/hour with XT Microslides\*



**The VITROS® XT 7600 Integrated System**

Up to 1,320 tests/hour with XT Microslide\*




**The VITROS® Automation Solution**

\*The numbers shown are maximum theoretical throughput. Actual throughput will vary by test mix and sample workflow.

Products not currently for sale in all countries

# VITROS® XT Systems delivers the consistently fastest results you can depend on.


## Fastest turnaround time for panel of common assays\*



\*TAT to complete the panel of these tests: Glucose, Calcium, Albumin, Total Protein, Sodium, Potassium, CO2, Chloride, BUN, Creatinine, ALP, ALT, AST, Bilirubin.

- VITROS® XT Solutions has an **up to 25% faster turnaround time** for common panel of tests compared to competitive systems.

## Results processed without requiring operator interventions



- VITROS® XT Solutions **reduce delays on 12.5% of samples** due to reviews/reruns/redraws compared to competitive systems.
- For labs doing 1MM tests per year, this provides **up to 300 more tests per day** without delay.

 **Consistently fastest. VITROS®.**

# Ensures your lab achieves the most accurate results with optimal system performance.

## Enabling accurate first-time results.

- **Accuracy:** minimize the impact of sample interferences with MicroSlides that filter out proteins and lipids
- **Confidence:** eliminate the risk of sample carryover and cross contamination with use of disposable tips
- **Trusted:** engineered to eliminate manual pre-screening for sample integrity with the MicroSensor's automatic hemolysis, icterus and turbidity testing
- **Improved:** results accuracy impacted by water quality variation

## Reliability your lab can depend on.

- **Trusted:** monitor, verify and document diagnostic checks during sample and assay processing with INTELLICHECK®
- **Improved:** fulfill more requests from difficult draws with sample volume as low as 2.7µl per test
- **Uptime:** 98% guaranteed uptime\*
- **Service:** fastest response time, same day service and remote resolution
- **Support:** e-Connectivity® minimized unplanned down time with secure, real-time, remote predictive monitoring and full engineer support

---

 **More accurate.**  **Best reliability. VITROS®.**

# More than a superior product. A comprehensive solution you can rely on.

## Ortho Care™ keeps your lab running smoothly

- **89%** same-day repair rate and **67%** remote resolution rate\*
- Predictive Technologies - World class web-based system monitoring software

**#1** Ranking by Service Trak™ for **superior service and support** received by Ortho Care™ - 4 years and counting

**#1** in **Overall Customer Satisfaction** for Clinical Chemistry, Immunoassay and Integrated Systems



**Successful ownership. VITROS®.**



\*Data on file.

# VITROS® XT Systems clears the path to deliver exceptional patient care.

## Lab results drive over 70% of clinical decisions.

When life is on the line, VITROS® XT is here for you with advanced technologies that propel your lab forward. And that makes all the difference in patient care.

## VITROS® XT delivers:

- Efficiency: exceptional productivity
- Accuracy: sophisticated interference-reduction technologies
- Reliability: outstanding 98% uptime guaranteed
- Ownership: integrated operation and service commitment

**...because every test is a life.**



[www.orthoclinicaldiagnostics.com](http://www.orthoclinicaldiagnostics.com)

### **About Ortho Clinical Diagnostics**

Ortho Clinical Diagnostics is a global leader of in vitro diagnostics serving the global clinical laboratory and immunohematology communities. Across hospitals, hospital networks, blood banks, and labs in more than 125 countries and territories, Ortho's high-quality products and services enable health care professionals to make better-informed treatment decisions. Ortho brings sophisticated testing technologies, automation, and information management and interpretation tools to clinical laboratories around the world to help them run more efficiently and effectively and improve patient care. For the immunohematology community, Ortho's blood typing products help ensure every patient receives blood that is safe, the right type and the right unit. Ortho's purpose is to improve and save lives with diagnostics, and it does that by reimagining what's possible. This is what has defined Ortho for more than 75 years, and it's what drives Ortho forward.

# Ortho Clinical Diagnostics